


Anti-Cdk6 Rabbit Monoclonal Antibody For Wb And Elisa - Q00534

Item Number: CM0002505



Introduction

High-quality rabbit monoclonal antibody against human CDK6. Validated for WB and ELISA, reacts with human and rat. Ideal for cell cycle and differentiation research. 37kDa.

[Learn More](#)

Parameter	Value
Product Name	CDK6 Rabbit mAb
Remarks/Alias	MCPH12; PLSTIRE; CDK6
Species	Human
GeneID (Human)	1021
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 227-326 of human CDK6 (Q00534).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q00534
Protein Weight	37kDa
Shipping	Ice bag