

Anti-Cdkn1A/P21Cip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P38936

Item Number: CM0024830



Introduction

Rabbit polyclonal antibody targeting human CDKN1A/p21CIP1 (p21). Validated for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse. Detects a band at ~18 kDa. Ideal for cancer and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	CDKN1A/p21CIP1 Rabbit pAb
Remarks/Alias	P21; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1; CDKN1A/p21CIP1
Species	Human
Gene ID (Human)	1026
Gene ID	1026
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human CDKN1A/p21CIP1 (NP_000380.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P38936
Protein Weight	18kDa
Shipping	Ice bag