

Anti-Cdx4 Polyclonal Antibody For Wb, Elisa - O14627

Item Number: CM0026259

Introduction

Rabbit polyclonal antibody against human CDX4, validated for WB and ELISA, cross-reacts with mouse. Ideal for studying homeobox protein function in developmental transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	CDX4 Rabbit pAb
Remarks/Alias	CDX4
Species	Human
Gene ID (Human)	1046
Gene ID	1046
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human CDX4 (NP_005184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14627
Protein Weight	30kDa
Shipping	Ice bag