


Anti-Cebpe Polyclonal Antibody For Wb, Elisa - Q15744

Item Number: CM0029699



Introduction

CEBPE rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human and mouse CEBPE (31kDa), a nuclear transcription factor critical for promyelocyte-myelocyte transition in myeloid differentiation.

[Learn More](#)

Parameter	Value
Product Name	CEBPE Rabbit pAb
Remarks/Alias	CRP1; SGD1; IMD108; C/EBP-epsilon; c/EBP epsilon; CEBPE
Species	Human
Gene ID (Human)	1053
Gene ID	1053
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human CEBPE (NP_001796.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15744
Protein Weight	31kDa
Shipping	Ice bag