

Anti-Cel/Bile Salt-Activated Lipase Polyclonal Antibody For Wb, Elisa - P19835

Item Number: CM0015488



Introduction

Rabbit polyclonal antibody against human CEL (bile salt-activated lipase), suitable for Western blot and ELISA with cross-reactivity to mouse. Key enzyme in dietary lipid absorption and fat-soluble vitamin metabolism.

[Learn More](#)

Parameter	Value
Product Name	CEL Rabbit pAb
Remarks/Alias	BAL; FAP; BSDL; BSSL; CELL; FAPP; LIPA; CEase; MODY8; CEL
Species	Human
Gene ID (Human)	1056
Gene ID	1056
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 350-550 of human CEL (NP_001798.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P19835
Protein Weight	79kDa
Shipping	Ice bag