


Anti-Celf1 Monoclonal Antibody For Wb, Ip, Elisa - Q92879

Item Number: CM0003926



Introduction

Rabbit monoclonal antibody against human CELF1 (CUGBP Elav-like family member 1). Suitable for WB, IP, ELISA. Cross-reacts with human, mouse, rat. Ideal for post-transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CELF1 Rabbit mAb
Remarks/Alias	CUGBP; NAB50; NAPOR; CUG-BP; CUGBP1; hNab50; BRUNOL2; EDEN-BP; CELF1
Species	Human
Gene ID (Human)	10658
Gene ID	10658
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CELF1 (Q92879).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92879
Protein Weight	52kDa
Shipping	Ice bag