


Anti-Cenpc Polyclonal Antibody For Wb And Elisa - Q03188

Item Number: CM0025660



Introduction

CENPC Polyclonal Antibody validated for WB and ELISA. Reacts with Human, Mouse, Rat. Ideal for research on kinetochore assembly, mitotic progression, and chromosome segregation.

[Learn More](#)

Parameter	Value
Product Name	CENPC Rabbit pAb
Remarks/Alias	MIF2, hcp-4, CENP-C, CENPC1, CENPC
Species	Human
Gene ID (Human)	1060
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CENPC (NP_001803.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03188
Protein Weight	107kDa
Shipping	Ice bag