

Anti-Cetn1 Rabbit Polyclonal Antibody For Wb, Elisa - Q12798

Item Number: CM0030812

Introduction

High-quality rabbit polyclonal antibody recognizing centrin-1 (CETN1) in human, mouse, and rat. Validated for Western blotting (WB) and ELISA; ideal for centrosome and microtubule-organizing center studies.

[Learn More](#)

Parameter	Value
Product Name	CETN1 Rabbit pAb
Remarks/Alias	CEN1; CETN; CETN1
Species	Human
Gene ID (Human)	1068
Gene ID	1068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CETN1 (NP_004057.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12798
Protein Weight	20kDa
Shipping	Ice bag