

Anti-Cgrrf1 Polyclonal Antibody For Wb, Elisa - Q99675

Item Number: CM0019004



Introduction

Rabbit polyclonal antibody targeting human CGRRF1 (CGR19/RNF197), validated for Western blot (WB) and ELISA applications. Suitable for cell growth regulation and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CGRRF1 Rabbit pAb
Remarks/Alias	CGR19; RNF197; CGRRF1
Species	Human
Gene ID (Human)	10668
Gene ID	10668
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-170 of human CGRRF1 (NP_006559.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99675
Protein Weight	38kDa
Shipping	Ice bag