

Anti-Chd3 Rabbit Polyclonal Antibody For Wb, Elisa - Q12873

Item Number: CM0023413

Introduction

Rabbit polyclonal antibody against CHD3 (Mi-2alpha), validated for WB and ELISA. For chromatin remodeling and NuRD complex research. Cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	CHD3 Rabbit pAb
Remarks/Alias	ZFH; Mi-2a; SNIBCPs; Mi2-ALPHA; CHD3
Species	Human
GeneID (Human)	1107
GeneID	1107
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1900-2000 of human CHD3 (NP_001005273.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12873
Protein Weight	227kDa
Shipping	Ice bag