


Anti-Chd4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14839

Item Number: CM0016346



Introduction

Rabbit polyclonal CHD4 antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. For chromatin remodeling and DNA repair studies. Swiss-Prot: Q14839.

[Learn More](#)

Parameter	Value
Product Name	CHD4 Rabbit pAb
Remarks/Alias	CHD-4; Mi-2b; SIHIWES; Mi2-BETA; CHD4
Species	Human
GeneID (Human)	1108
GeneID	1108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1520-1690 of human CHD4 (NP_001264.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14839
Protein Weight	218kDa
Shipping	Ice bag