

Anti-Chk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O96017

Item Number: CM0016467



Introduction

Rabbit polyclonal antibody against human Chk2 (CHEK2) kinase, validated for WB, IF/ICC, IP, and ELISA. Ideal for studying DNA damage response, cell cycle checkpoints, and apoptosis signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	Chk2 Rabbit pAb
Remarks/Alias	CDS1; CHK2; LFS2; RAD53; hCds1; HuCds1; PP1425
Species	Human
Gene ID (Human)	11200
Gene ID	11200
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human Chk2 (NP_009125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	O96017
Protein Weight	61kDa
Shipping	Ice bag