


Anti-Chmp6 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Fz7

Item Number: CM0029974



Introduction

CHMP6 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA; detects human and mouse CHMP6, an ESCRT-III component involved in endosomal sorting and MVB formation. Ideal for research on membrane trafficking and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	CHMP6 Rabbit pAb
Remarks/Alias	VP520; CHMP6
Species	Human
Gene ID (Human)	79643
Gene ID	79643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-200 of human CHMP6 (NP_078867.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96FZ7
Protein Weight	23kDa
Shipping	Ice bag