

Anti-Chp1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99653

Item Number: CM0020595



Introduction

Rabbit polyclonal antibody against CHP1, validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets calcineurin B homologous protein 1 (22 kDa), involved in pH regulation and vesicular trafficking.

[Learn More](#)

Parameter	Value
Product Name	CHP1 Rabbit pAb
Remarks/Alias	CHP; p22; p24; SPAX9; Sid470p; SLC9A1BP; CHP1
Species	Human
GeneID (human)	11261
GeneID	11261
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human CHP1 (NP_009167.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99653
Protein Weight	22kDa
Shipping	Ice bag