

Anti-Chrng Polyclonal Antibody For Wb, Elisa - P07510

Item Number: CM0028088

Introduction

High-quality rabbit polyclonal antibody targeting the human acetylcholine receptor subunit gamma (CHRNA7). Validated for Western blot and ELISA, with cross-reactivity to human, mouse, and rat samples. Suitable for ion channel and neuromuscular junction research applications.

[Learn More](#)

Parameter	Value
Product Name	CHRNA7 Rabbit pAb
Remarks/Alias	ACHR7; CHRNA7
Species	Human
GeneID (Human)	1146
GeneID	1146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 328-517 of human CHRNA7 (NP_005190.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07510
Protein Weight	58kDa
Shipping	Ice bag