


Anti-Ciao1 Rabbit Polyclonal Antibody For Wb, Elisa - 076071

Item Number: CM0030507



Introduction

CIAO1 rabbit polyclonal antibody validated for WB and ELISA, targeting human CIAO1 (CIA1/WDR39), a cytosolic iron-sulfur assembly factor; cross-reacts with mouse and rat. Ideal for Fe/S protein biogenesis and chromosome segregation research.

[Learn More](#)

Parameter	Value
Product Name	CIAO1 Rabbit pAb
Remarks/Alias	CIA1; WDR39; CIAO1
Species	Human
Gene ID (Human)	9391
Gene ID	9391
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human CIAO1 (NP_004795.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076071
Protein Weight	38kDa
Shipping	Ice bag