


Anti-Ciap2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13489

Item Number: CM0017133



Introduction

Rabbit polyclonal antibody targeting human cIAP2 (BIRC3, API2), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for apoptosis, NF-kB signaling, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	cIAP2 Rabbit pAb
Remarks/Alias	AIP1; API2; MIHC; CIAP2; HAIP1; HIAP1; IAP-1; MALT2; RNF49; c-IAP2; cIAP2
Species	Human
Gene ID (Human)	330
Gene ID	330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-130 of human cIAP2 (NP_001156.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13489
Protein Weight	68kDa
Shipping	Ice bag