


Anti-Cib1 Polyclonal Antibody For Wb, Elisa - Q99828

Item Number: CM0028384



Introduction

Rabbit polyclonal antibody targeting human CIB1, validated for Western blotting (WB) and ELISA. Reacts with human, mouse, rat samples. Suitable for studying cell adhesion, integrin signaling, and cancer biology.

[Learn More](#)

Parameter	Details
Product Name	CIB1 Rabbit pAb
Remarks/Alias	CIB; CIBP; KIP1; PRKDCIP; SIP2-28; CIB1
Species	Human
GeneID (Human)	10519
GeneID	10519
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human CIB1 (NP_006375.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99828
Protein Weight	22kDa
Shipping	Ice bag