

Anti-Cip4 Rabbit Monoclonal Antibody For Wb And Elisa - Q15642

Item Number: CM0002118



Introduction

Cip4 rabbit monoclonal antibody validated for WB and ELISA. Reacts with human and mouse Cdc42-interacting protein 4, involved in endocytosis and insulin signaling. Suitable for cytoskeleton and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	Cip4 Rabbit mAb
Remarks/Alias	STP; CIP4; HSTP; STOT; TRIP-10; Cip4
Species	Human
GeneID	9322
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cip4 (Q15642).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15642
Protein Weight	68kDa
Shipping	Ice bag