


## Anti-Ckap5 Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, If-P, Elisa - Q14008

Item Number: CM0024828



### Introduction

CKAP5 rabbit polyclonal antibody raised against human CKAP5 (ch-TOG), validated for WB, IP, IF/ICC, IF-P, and ELISA with cross-reactivity to mouse and rat. Ideal for investigating microtubule dynamics, mitotic spindle organization, and tumor biology.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CKAP5 Rabbit pAb                                                                                                         |
| Remarks/Alias    | TOG, MSPS, TOGp, CHTOG, ch-TOG, CKAP5                                                                                    |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 9793                                                                                                                     |
| GeneID           | 9793                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1763-2032 of human CKAP5 (NP_001008938.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IP, IF/ICC, IF-P, ELISA                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | Q14008                                                                                                                   |
| Protein Weight   | 225kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |