

## Anti-Clathrin Heavy Chain Rabbit Monoclonal Antibody For Wb, Ihc, If, Ip, Elisa - Q00610 / P53675

Item Number: CM0010001



### Introduction

High-specificity rabbit monoclonal antibody against human Clathrin heavy chain (CLTC/CLTCL1; Q00610/P53675). Validated for WB, IHC, IF, IP, ELISA. Reacts with human, mouse, and rat. Ideal for endocytosis, mitosis, and autophagy research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Clathrin heavy chain Rabbit mAb                                                                                                       |
| Remarks/Alias    | Hc; CHC; CHC17; MRD56; CLH-17; CLTCL2; Clathrin heavy chain                                                                           |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 1213/8218                                                                                                                             |
| Gene ID          | 1213/8218                                                                                                                             |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Clathrin heavy chain (Q00610). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | Q00610/P53675                                                                                                                         |
| Protein Weight   | 191kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                               |