

Anti-Claudin-3 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - O15551

Item Number: CM0026475

Introduction

High-quality rabbit polyclonal antibody against human Claudin-3 (CLDN3). Validated for Western blot, immunofluorescence, and ELISA with cross-reactivity to mouse and rat. Ideal for tight junction research.

[Learn More](#)

Parameter	Value
Product Name	CLDN3 Rabbit pAb
Remarks/Alias	RVP1; HRVP1; C7orf1; CPE-R2; CPETR2; CLDN3
Species	Human
GenelD	1365
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 141-220 of human CLDN3 (NP_001297.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15551
Protein Weight	23kDa
Shipping	Ice bag