

Anti-Claudin18.2 Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P56856-2

Item Number: CM0011266

Introduction

High-quality anti-Claudin18.2 rabbit monoclonal antibody validated for IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for gastric tight junction and barrier function research. SWISS: P56856-2.

[Learn More](#)

Parameter	Specification
Product Name	Claudin18.2 Rabbit mAb
Remarks/Alias	SFTA5; SFTPJ
Species	Human
GenelD (Human)	51208
GenelD	51208
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56856-2
Protein Weight	28kDa
Shipping	Ice bag