


Anti-Clcnka Rabbit Polyclonal Antibody For Wb, Elisa - P51800

Item Number: CM0028531



Introduction

CLCNKA Rabbit Polyclonal Antibody validated for WB and ELISA in human and mouse samples. Targets the chloride channel CIC-Ka, involved in hearing and renal function.

[Learn More](#)

Parameter	Value
Product Name	CLCNKA Rabbit pAb
Remarks/Alias	CLCK1; CIC-K1; hCIC-Ka; CLCNKA
Species	Human
Gene ID (Human)	1187
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-644 of human CLCNKA (NP_001244068.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51800
Protein Weight	75kDa
Shipping	Ice bag