

## Anti-Clon11 Rabbit Polyclonal Antibody For Wb, Elisa - 075508

Item Number: CM0029507

### Introduction

High-quality CLDN11 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse claudin-11, a tight junction protein involved in cell adhesion and leukodystrophy.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | CLDN11 Rabbit pAb                                                                                                  |
| Remarks/Alias    | OSP; OTM; HLD22; CLDN11                                                                                            |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 5010                                                                                                               |
| GeneID           | 5010                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-207 of human CLDN11 (NP_005593.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | 075508                                                                                                             |
| Protein Weight   | 22kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |