

Anti-Clndnd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ny35

Item Number: CM0025632

Introduction

Rabbit polyclonal antibody targeting human CLDND1, a tight junction protein regulating small molecule permeability. Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for tight junction research.

[Learn More](#)

Parameter	Value
Product Name	CLDND1 Rabbit pAb
Remarks/Alias	Z38; C3orf4; GENX-3745; CLDND1
Species	Human
GeneID (Human)	56650
GeneID	56650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-138 of human CLDND1 (Q9NY35).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NY35
Protein Weight	29kDa
Shipping	Ice bag