


Anti-Cleaved Caspase-8 (Asp384) Rabbit Monoclonal Antibody For Wb, Elisa - Q14790

Item Number: CM0004729



Introduction

Rabbit monoclonal antibody targeting cleaved Caspase-8 (Asp384). Validated for Western blot and ELISA in human samples. Suitable for apoptosis, necroptosis, and pyroptosis research applications.

[Learn More](#)

Parameter	Value
Product Name	Cleaved Caspase-8 (Asp384) Rabbit mAb
Remarks/Alias	CAP4; MACH; MCH5; FLICE; ALPS2B; Casp-8
Species	Human
GeneID (Human)	841
GeneID	841
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 371-470 of human Caspase-8 (NP_203519.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14790
Protein Weight	55kDa
Shipping	Ice bag