


Anti-CLN5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075503

Item Number: CM0013159



Introduction

CLN5 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Recognizes bis(monoacylglycero)phosphate synthase CLN5, a lysosomal enzyme involved in BMP synthesis and neuronal ceroid lipofuscinosis. Ideal for neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	CLN5 Rabbit pAb
Remarks/Alias	CLN5; NCL
Species	Human
GeneID (human)	1203
GeneID	1203
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-407 of human CLN5 (NP_006484.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075503
Protein Weight	41kDa
Shipping	Ice bag