

Anti-Clns1A Polyclonal Antibody For Wb, Ihc-P, Elisa - P54105

Item Number: CM0019517



Introduction

Rabbit polyclonal antibody against human CLNS1A (pICln) for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Involved in spliceosomal snRNP assembly and Sm protein methylation. Suitable for RNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	CLNS1A Rabbit pAb
Remarks/Alias	CLCI, ICln, CLNS1B, CLNS1A
Species	Human
GeneID (Human)	1207
GeneID	1207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human CLNS1A (NP_001284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54105
Protein Weight	26kDa
Shipping	Ice bag