

Anti-CLNS1A Rabbit Polyclonal Antibody For Wb, Elisa - P54105

Item Number: CM0024834



Introduction

Rabbit polyclonal antibody targeting human CLNS1A (methylosome subunit pICln) for use in Western blot and ELISA. Validated in human samples; essential for spliceosome assembly and mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	CLNS1A Rabbit pAb
Remarks/Alias	CLCI; ICln; CLNS1B; CLNS1A
Species	Human
GeneID (Human)	1207
GeneID	1207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human CLNS1A (NP_001284.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P54105
Protein Weight	26kDa
Shipping	Ice bag