

Anti-Clptm1 Rabbit Monoclonal Antibody For Wb And Elisa - O96005

Item Number: CM0011077

Introduction

Rabbit monoclonal antibody against human CLPTM1 (cleft lip and palate transmembrane protein 1). Suitable for Western blot and ELISA, cross-reacts with human. Research use for GPI-anchored protein processing and GABAergic signaling.

[Learn More](#)

Parameter	Value
Product Name	CLPTM1 Rabbit mAb
Remarks/Alias	CLPTM1; transmembrane protein
Species	Human
Gene ID (Human)	1209
Gene ID	1209
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CLPTM1 (O96005)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O96005
Protein Weight	76kDa
Shipping	Ice bag