


Anti-Clspn (Claspin) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Haw4

Item Number: CM0023112



Introduction

Rabbit polyclonal antibody targeting human CLSPN/Claspin (UniProt Q9HAW4). Suitable for WB and ELISA. Ideal for investigating DNA damage checkpoints, replication fork stability, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	CLSPN Rabbit pAb
Remarks/Alias	CLSPN; claspin
Species	Human
GeneID (human)	63967
GeneID	63967
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 650-750 of human CLSPN (NP_001177410.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HAW4
Protein Weight	151kDa
Shipping	Ice bag