

Anti-Cog2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14746

Item Number: CM0027737

Introduction

CamelBio offers a rabbit polyclonal antibody targeting COG2 (Q14746), a conserved oligomeric Golgi complex subunit essential for Golgi morphology and function. Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for Golgi trafficking and congenital disorders of glycosylation research.

[Learn More](#)

Parameter	Value
Product Name	COG2 Rabbit pAb
Remarks/Alias	LDLC; CDG2Q; COG2
Species	Human
GeneID (Human)	22796
GeneID	22796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human COG2 (NP_031383.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14746
Protein Weight	83kDa
Shipping	Ice bag