

Anti-Cog3 Polyclonal Antibody For Wb, Elisa - Q96Jb2

Item Number: CM0021386

Introduction

Rabbit polyclonal antibody targeting human COG3 (SEC34), validated for Western blot (WB) and ELISA. Reacts with human, mouse, rat. Useful for Golgi transport studies. 94 kDa.

[Learn More](#)

Parameter	Value
Product Name	COG3 Rabbit pAb
Remarks/Alias	SEC34; COG3
Species	Human
Gene ID	83548 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human COG3 (NP_113619.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96JB2
Protein Weight	94kDa
Shipping	Ice bag