

Anti-Complement C3 Rabbit Polyclonal Antibody For If-P, Elisa - P01024

Item Number: CM0013238



Introduction

Rabbit polyclonal antibody against human Complement C3 (aa 1364-1663). Validated for IF-P and ELISA; cross-reacts with human, mouse, and rat. Suitable for complement pathway and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Complement C3 Rabbit pAb
Remarks/Alias	ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p; Complement C3
Species	Human
Gene ID (Human)	718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1364-1663 of human Complement C3 (NP_000055.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01024
Protein Weight	187kDa
Shipping	Ice bag