


Anti-Connexin 43 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P17302

Item Number: CM0013551



Introduction

High-quality rabbit polyclonal antibody targeting human Connexin 43 (Cx43, GJA1). Validated for WB, IF/ICC, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for gap junction, cardiac, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Connexin 43 Rabbit pAb
Remarks/Alias	HSS; CMDR; CX43; EKVP; GJAL; ODDD; AVSD3; EKVP3; HLHS1; PPKCA; Connexin 43
Species	Human
Gene ID (Human)	2697
Gene ID	2697
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 233-382 of human Connexin 43 (NP_000156.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17302
Protein Weight	43 kDa
Shipping	Ice bag