


Anti-Copa Polyclonal Antibody For Wb And Elisa - P53621

Item Number: CM0021536



Introduction

Rabbit polyclonal antibody against human COPA (coatamer subunit alpha). Validated for Western blot and ELISA; cross-reacts with mouse. Ideal for studying COPI vesicle trafficking and Golgi function.

[Learn More](#)

Parameter	Specification
Product Name	COPA Rabbit pAb
Remarks/Alias	AILJK; HEP-COP; alpha-COP; COPA
Species	Human
GeneID	1314
Immunogen	Synthetic peptide (aa 550-650 of human COPA, NP_001091868.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P53621
Protein Weight	138kDa
Shipping	Ice bag