

Anti-Coq3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nzj6

Item Number: CM0029810

Introduction

Anti-COQ3 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Targets mitochondrial ubiquinone biosynthesis O-methyltransferase. Cross-reacts with human, mouse, and rat COQ3. Ideal for ubiquinone pathway research.

[Learn More](#)

Parameter	Value
Product Name	COQ3 Rabbit pAb
Remarks/Alias	DHHBMT; bA9819.1; DHHBMTASE; UG0215E05; COQ3
Species	Human
Gene ID (Human)	51805
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 50-369 of human COQ3 (NP_059117.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZJ6
Protein Weight	41kDa
Shipping	Ice bag