

Anti-Corin Polyclonal Antibody For Wb, Elisa - Q9Y5Q5

Item Number: CM0018035

Introduction

Rabbit polyclonal antibody against human CORIN, suitable for WB and ELISA. Recognizes human, mouse, and rat. Recombinant immunogen (aa 1-110). For cardiovascular and natriuretic peptide research.

[Learn More](#)

Parameter	Value
Product Name	CORIN Rabbit pAb
Remarks/Alias	CRN; ATC2; Lrp4; PEE5; TMPRSS10; CORIN
Species	Human
Gene ID (Human)	10699
Gene ID	10699
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human CORIN (NP_006578.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5Q5
Protein Weight	116kDa
Shipping	Ice bag