

Anti-Coro1A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31146

Item Number: CM0026223

Introduction

CORO1A Rabbit Polyclonal Antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat coronin-1A, a 51 kDa actin-binding protein involved in phagosome maturation and cell motility. Ideal for immunology and infection research.

[Learn More](#)

Parameter	Value
Product Name	CORO1A Rabbit pAb
Remarks/Alias	p57, IMD8, TACO, CLABP, HCORO1, CLIPINA, CORO1A
Species	Human
Gene ID (Human)	11151
Gene ID	11151
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 362-461 of human CORO1A (NP_001180262.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31146
Protein Weight	51kDa
Shipping	Ice bag