

Anti-Cortactin Polyclonal Antibody For Wb, If/Icc, Elisa - Q14247

Item Number: CM0029247



Introduction

Rabbit polyclonal antibody targeting human cortactin (EMS1), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for cytoskeleton and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	Cortactin Rabbit pAb
Remarks/Alias	EMS1; Cortactin
Species	Human
GeneID (human)	2017
GeneID	2017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-400 of human Cortactin (NP_005222.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14247
Protein Weight	62kDa
Shipping	Ice bag