


Anti-Cox Iv Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P19783

Item Number: CM0029255



Introduction

High-specificity rabbit polyclonal antibody targeting COX IV (Cytochrome c oxidase subunit 4 isoform 1) validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for mitochondrial respiration and oxidative phosphorylation studies.

[Learn More](#)

Parameter	Value
Product Name	COX IV Rabbit pAb
Remarks/Alias	COX4; COXIV; COX4-1; COXIV-1; MC4DN16; COX IV-1
Species	Mouse
Gene ID (Human)	1327
Gene ID	12857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-98 of mouse COX IV (NP_001280488.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19783
Protein Weight	20 kDa
Shipping	Ice bag