

## Anti-Cox7A1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P24310

Item Number: CM0003815

### Introduction

High-quality anti-COX7A1 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat COX7A1. Ideal for studying mitochondrial complex IV assembly and respiratory supercomplexes.

[Learn More](#)

| Parameter        | Value                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------|
| Product Name     | COX7A1 Rabbit mAb                                                                                      |
| Remarks/Alias    | COX7A; COX7AH; COX7AM; COX7A1                                                                          |
| Species          | Human                                                                                                  |
| GeneID (Human)   | 1346                                                                                                   |
| GeneID           | 1346                                                                                                   |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-79 of human COX7A1 (NP_001855.1). |
| Source           | Rabbit                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                  |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                      |
| SWISS            | P24310                                                                                                 |
| Protein Weight   | 9kDa                                                                                                   |
| Shipping         | Ice bag                                                                                                |