

Anti-Cox7A2L Rabbit Polyclonal Antibody For Wb, Elisa - O14548

Item Number: CM0027949

Introduction

COX7A2L rabbit polyclonal antibody, validated for WB and ELISA, detects human, mouse, and rat COX7A2L (SCAF1). Recombinant immunogen targeting N-terminal region. Ideal for studying mitochondrial respiratory supercomplex assembly and energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	COX7A2L Rabbit pAb
Remarks/Alias	EB1, SCAF1, SCAFI, SIG81, COX7AR, COX7RP, COX7A2L
Species	Human
GeneID (Human)	9167
GeneID	9167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human COX7A2L (NP_004709.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14548
Protein Weight	13kDa
Shipping	Ice bag