

Anti-Cox8A Rabbit Polyclonal Antibody For Wb And Elisa - P10176

Item Number: CM0025866



Introduction

COX8A Rabbit Polyclonal Antibody validated for WB, ELISA. Reacts with Human, Mouse. Targets Cytochrome c oxidase subunit 8A, mitochondrial (P10176), 8kDa. Ideal for mitochondrial respiratory chain studies.

[Learn More](#)

Parameter	Value
Product Name	COX8A Rabbit pAb
Remarks/Alias	COX; COX8; VIII; COX8L; COX8-2; VIII-L; MC4DN15; COX8A
Species	Human
GeneID (Human)	1351
GeneID	1351
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human COX8A (NP_004065.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10176
Protein Weight	8kDa
Shipping	Ice bag