

Anti-Cpa4 (Carboxypeptidase A4) Polyclonal Antibody For Wb, Elisa - Q9Ui42

Item Number: CM0021722



Introduction

Rabbit polyclonal antibody against human CPA4 (carboxypeptidase A4), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying metalloprotease function and peptide hormone catabolism. Research use only.

[Learn More](#)

Parameter	Value
Product Name	CPA4 Rabbit pAb
Remarks/Alias	CPA3; CPA4
Species	Human
GeneID (human)	51200
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 114-280 of human CPA4 (NP_057436.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Ui42
Protein Weight	47kDa
Shipping	Ice bag