

Anti-Cpe Rabbit Polyclonal Antibody For Wb And Elisa - P16870

Item Number: CM0028238

Introduction

High-quality rabbit polyclonal antibody against carboxypeptidase E (CPE). Validated for WB and ELISA, cross-reacts with human and mouse. Supports research in neuroendocrine biology and hormone processing.

[Learn More](#)

Parameter	Details
Product Name	CPE Rabbit pAb
Remarks/Alias	CPH; BDVS; IDDHH; CPE
Species	Human
Gene ID	1363
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 257-476 of human CPE (NP_001864.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16870
Protein Weight	53kDa
Shipping	Ice bag