

Anti-Cr2/Cd21 Monoclonal Antibody For Wb, Ihc-P, Elisa - P20023

Item Number: CM0006765



Introduction

Rabbit monoclonal antibody targeting human CR2/CD21 (Complement receptor type 2), validated in WB, IHC-P, ELISA. Ideal for studying B-cell activation, complement pathway, and EBV receptor interactions.

[Learn More](#)

Parameter	Value
Product Name	CR2/CD21 Rabbit mAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9; CR2/CD21
Species	Human
GeneID (Human)	1380
GeneID	1380
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 934-1033 of human CR2/CD21 (NP_001868.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P20023
Protein Weight	113kDa
Shipping	Ice bag