

Anti-Cr2/Cd21 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P20023

Item Number: CM0007692



Introduction

Rabbit monoclonal antibody against human CR2/CD21, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for B-cell receptor, complement pathway, and Epstein-Barr virus research.

[Learn More](#)

Parameter	Value
Product Name	CR2/CD21 Rabbit mAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9; CR2/CD21
Species	Human
Gene ID (Human)	1380
Gene ID	1380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 933-1033 of human CR2/CD21. (NP_001868.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20023
Protein Weight	113kDa
Shipping	Ice bag