

Anti-Creb3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43889

Item Number: CM0026118

Introduction

Rabbit polyclonal antibody against human CREB3 (LZIP/LUMAN), validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat CREB3. Targets ER stress transcription factor involved in UPR, immunity, and viral pathogenesis.

[Learn More](#)

Parameter	Value
Product Name	CREB3 Rabbit pAb
Remarks / Alias	LZIP; LUMAN; sLZIP; CREB3
Species	Human
GeneID (Human)	10488
GeneID	10488
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human CREB3 (NP_006359.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43889
Protein Weight	41kDa
Shipping	Ice bag